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
This manual should be available for all users of these equipments. To get the best results and a higher duration of this equipment it is advisable to read carefully this manual and follow the processes of use.

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AUXILAB S.L.

**DECLARATION OF CONFORMITY
NAHITA MICROTOME of Auxilab, S.L for the Directive
of Machines (89/392/CEE modified) and the regulations
adopted for their transposition**

NAME OF THE MANUFACTURER / IMPORTER:

AUXILAB, S.L.

ADDRESS:

**Polígono Morea Norte, 8
31191 Beriain (Navarra)**

WE STATE THAT:

**MICROTOME 591
CODE 50591000**

Is designed and manufactured according to:

- ◆ Directive 89/392/CEE, including the modifications and national regulations that transpose them.
- ◆ And that the following harmonized rules have been applied (or a part of them):
UNE 61010-1/A2, UNE-EN61010-2-051.

BERIAIN 28th February 2007

Signed by: ALFONSO AINCIBURU SANZ
DIRECTOR/MANAGER

Polígono Morea Norte, 8 31191 Beriain (Navarra) - Spain. Tel. 948 310 513 Fax 948 312 071
Internet: www.auxilab.es · Email: correo@auxilab.es

TROUBLE	CAUSE	SOLUCIÓN
	Knife clearance angle too wide.	Systematically reduce clearance angle until optimum angle width has been found.
	Insufficient clamping of sample and/or knife holder	Check all screws and clamping mechanisms at sample holder system and knife holder. Retighten if necessary.
No blade feed, no more sections	Blade has reached front limit stop position.	Move blade backwards via coarse feed knob.
Blades wear out very quickly	Cutting forces too high	Select sectioning speed or section thickness. Reduce section thickness or turn handwheel at slower speed.

INSTRUCTIONS ON ENVIRONMENT PROTECTION

At the end of its life cycle, please, do not dispose of this equipment by throwing it in the usual garbage; hand it over a collection point for the recycling of electrical and electronic appliances. It does not contain dangerous or toxic products for humans but a non adequate disposal would damage the environment.

The materials are recyclable as mentioned in its marking. By recycling material or by other forms of re-utilization of old appliances, you are making an important contribution to protect our environment.

Please inquire at the community administration for the authorized disposal location.

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
Cleaning

- ◆ Never use scourers or substances that can grate for cleaning metallic parts such as stainless steel, aluminium, coatings, etc. as they damage the microtome and produce an early ageing of the equipment. Neither use alcohol or alcohol containing detergents nor abrasive detergents or acetone or xylene containing solvents.
- ◆ Use a fluff-free cloth dampened with soaped water that does not contain abrasives.

ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.

7. TROUBLESHOOTING

TROUBLE	CAUSE	SOLUCIÓN
Alternately thick and thin sections are produced In extreme cases, no section is obtained.	Blade insufficiently clamped	Re-clamp blade
	Blunt cutting edge	Adjust knife holder laterally or use new blade.
	Knife clearance angle too narrow.	Systematically try wider clearance angle until optimum angle width has been found.
Sections are extremely compressed, wrinkled or jammed together	Blunt cutting edge	Use different zone of cutting blade or use new blade.
	Sample too warm	Cool sample prior to sectioning
	Sectioning speed too high	Work at lower sectioning speed
Noise during sectioning Knife "sing." when sectioning hard samples. Sections show scratches.	Sectioning speed too high	Work at lower sectioning speed

- ◆ Put the microtome on top of a horizontal, plane and stable table making a free space at least at 30 cm per side.
- ◆ Do not expose the microtome directly to the sun, etc. Avoid vibrations, dust and dry environments.
- ◆ Always remove the blade before detaching the knife holder from the instrument.
- ◆ Never try to catch a falling knife.
- ◆ Never place a knife anywhere with the cutting edge facing upwards.
- ◆ Always clamp the sample block before clamping the knife.
- ◆ When working with cleaning detergents, please bear in mind manufacturer safety instructions as well as current laboratory regulations.
- ◆ Wear glasses in case of working with samples that can produce splinters.
- ◆ Always lock the handwheel prior to equipment cleaning (see Cleaning), prior to any manipulation of knife or sample and during all work breaks.
- ◆ No fluid must enter the interior of the equipment during cleaning.
- ◆ Never try to repair the microtome by yourself, since you will lose the warranty and may provoke damages to the general operating system.
- ◆ Made under the European regulations for security on machines.

6. MAINTENANCE AND CLEANING


To get the best results and a higher duration of the microtome it is essential to follow the processes of use.

Note: All the processes of use mentioned below will not have any value unless you keep a continued and careful maintenance.

- ◆ Please follow the processes of use of this manual.
- ◆ This manual should be available for all users of this equipment.
- ◆ Always use original components and supplies. Other devices can be similar but they can damage the equipment.
- ◆ Never try to repair the microtome by yourself, since you will lose the warranty and may provoke damages to the general operating system.
- ◆ In the event of breakdown please contact your distributor to overhaul through Nahita Technical Assistance Department.

Before each cleaning:

- ◆ Move sample clamp to uppermost position and lock handwheel.
- ◆ Insert the knife into the receptacle for used blades in the base of the blade dispenser.
- ◆ Remove knife holder and knife holder base for cleaning.
- ◆ Remove sample from sample clamp.
- ◆ Remove section waste with a dry brush.
- ◆ Remove sample clamp and clean it.

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5. MODO DE USO

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- ◆ Orientate sample in north-south direction with screw (1.16) and in east-west direction with screw (1.17).
- ◆ Turn clamp upwards.

5.5 Trimming the sample

Sample can be trimmed via the section thickness selection (max. 25 µm):

- ◆ Approach the knife holder to the sample via the coarse feed knob (1.8).
- ◆ Select section thickness (e.g. 25 µm) with the adjusting knob (1.7).
- ◆ Turn handwheel (1.10) until sample has been trimmed down to the desired level.

5.6 Sectioning

NOTE: always rotate handwheel evenly. Rotating speed must be adjusted to the sample hardness; harder samples should be cut at slower speed.

Never cut sample with the same region of the blade edge used for trimming. The knife holder base presents a lateral displacement system that allows the use of the entire width of the blade without having to release the knife:

- ◆ Turn backwards the clamp located at the left side of the knife holder base (1.2) to release.

- ◆ Move laterally the knife holder base with the knife holder to the desired position.

Hold the knife holder base at this position by turning forwards the lever (1.2).

Before sectioning it is essential to unlock the manual handwheel; for this purpose pull the locking device and give a half-turn to the handwheel to place it on cutting position.

- ◆ Select the section thickness (in microns) through the corresponding knob (1.7)

For sectioning, turn handwheel evenly (1.10) in clockwise direction.

IMPORTANT: Do not lock the coarse feed knob while rotating the manual handwheel. If the coarse feed knob is locked the knife does not advance and sectioning is not possible.

◆

Pick up the obtained sections and mount them on microscope slides.

5.7 Finishing operation

- ◆ Move sample clamp to uppermost position and lock the handwheel.
 - ◆ Insert the knife into the receptacle for used blades in the base of the blade dispenser.
 - ◆ Remove sample from sample clamp.
 - ◆ Remove all section waste.
- Clean the equipment (see Cleaning).

5.8 Security

The microtome must be used by previously qualified staff that knows how the equip-

wheel must be always locked before changing the sample as well as during all work breaks.

- ◆ Move the holding device for sample clamp (1.12) to the uppermost position; for this purpose, turn the manual handwheel (1.10) and lock it by the device located at the lower side (1.11).
- ◆ Insert the cassette horizontally or vertically into the sample clamp. For this purpose, pull forward the expelling lever (1.15) and place the cassette. Release the lever to secure the cassette in position.

5.2 Placing the knife

Nahita microtome 591 is equipped with a knife holder for standard disposable blades of common use.

NOTE: *Be careful when handling disposable blades for microtome. The blade edge is extremely sharp and can cause severe injuries*

- ◆ Turn upwards the knife holding clamp (1.5).
- ◆ Carefully insert the disposable blade from one side into the knife holder.
- ◆ To hold the knife, turn downwards the holding clamp (1.5).
- ◆ Pay attention to the blade edge; replace the knife when the obtained sections present scratches that could make them useless.

5.3 Adjusting the clearance angle

The marks (0°, 5° and 10°) for clearance angle adjustment (1.6) are on the right side of the knife holder. Another mark is on the right side of the knife holder base which serves as a reference point to select the desired angle.

- ◆ Rotate anticlockwise the clamping lever for knife holder (1.3) to release.
- ◆ Adjust the knife holder position until the selected mark matches the reference mark on the knife holder base.
- ◆ Hold the knife holder at this position and turn clockwise the clamping lever (1.3) to secure it.

5.4 Sample orientation

Orientation system of the sample holder allows adjusting the position of the block surface without having to release the clamping lever.

- ◆ Unlock the coarse feed knob (1.8) by turning to the right the locking lever (1.9).
- ◆ Rotate the coarse feed knob and make advance the knife holder base with the knife holder until nearly reaching the sample.
- ◆ Move the holding device for sample clamp (1.12) to the uppermost position and lock the manual handwheel.
- ◆ Release the locking lever of the orientation system (1.18) by turning downwards.

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3. TECHNICAL SPECIFICATIONS

<i>Slice thickness</i>	1-25 µm
<i>Min. adjustment</i>	1 µm
<i>Clearance angle</i>	0-10°
<i>Coarse feed</i>	Manual, through coarse feed knob
<i>Trimming thickness</i>	Adjustable max. 25 µm
<i>Dimensions (LxWxH)</i>	320x380x320 mm
<i>Weight</i>	42 Kg (approx.)

4. SETTING-UP

Unwrap the microtome and take off all the protections in which it comes fitted. Make sure that the equipment does not present any damage because of the shipment. In case it presents any damage, tell it immediately to your transport agent or dealer so that they can make the claims in the correct time limit.

Please keep the original wrapping; you will always need it for returns enclosed with all the accessories supplied.

We will only accept any equipment return within 15 days after delivery and provided it comes in its original wrapping.

4.1 Mounting the universal cassette clamp

The holding device of the sample holder will let you adjust the cassette clamp:

- ◆ Move the holding device for sample clamp (1.12) to the uppermost position; for this purpose, turn the manual handwheel (1.10) and lock it by the device located at the lower side (1.11).
- ◆ For placing the cassette clamp turn anticlockwise the clamping lever for sample clamp (1.13).
- ◆ Insert the cassette clamp (1.14) from the right side and slide it.
- ◆ To secure the cassette clamp, turn clockwise the clamping lever (1.13).

4.2 Mounting the knife holder


- ◆ Turn anticlockwise the camping lever for knife holder (1.3) to release it.
- ◆ Slide the groove of knife holder (1.4) onto T-piece of knife holder base.
- ◆ Hold the knife holder by turning clockwise the corresponding camping lever (1.3).

5. OPERATION

5.1 Placing the sample

NOTE: Always clamp the sample block before clamping the knife. The manual hand

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Knife holder: the approach of the knife towards the sample can be adjusted through a coarse feed knob.

Advance mechanism: in rotary microtomes the sample holder advances over the blade through an up-down movement. A complex mechanical system is necessary to perform these movements that finally allow obtaining serial tissue sections from the paraffin block.

Knife: universal disposable blades.

Manual handwheel: it transmits to the advance mechanism the impulse necessary for its correct operation; each complete turn of the handwheel produces a cut sequence. The handwheel presents a safety device that blocks it for its transport or when it is not being used. So, before using the microtome it is essential to pull the blocking mechanism and give a half-turn to the handwheel to place it on cutting position.

- | | |
|--|---|
| 1. Knife holder base | 10. Manual handwheel |
| 2. Clamping lever for lateral displacement | 11. Handwheel locking device |
| 3. Clamping lever for knife holder | 12. Holding device for sample clamp |
| 4. Knife holder | 13. Clamping lever for sample clamp |
| 5. Clamping lever for knife | 14. Universal cassette clamp |
| 6. Clearance angle scale | 15. Cassette expelling lever |
| 7. Slice thickness adjusting knob | 16. North/South orientation screw |
| 8. Coarse feed knob | 17. East/West orientation screw |
| 9. Coarse feed locking lever | 18. Locking lever of the orientation system |

INDEX OF LANGUAGES

Spanish	2-12
English.....	13-23

INDEX OF CONTENTS

1. USES OF THE INSTRUMENT	14
2. DESCRIPTION.....	14
3. TECHNICAL SPECIFICATIONS	17
4. SETTING UP	17
5. OPERATION	17
6. MAINTENANCE AND CLEANING	20
7. TROUBLESHOOTING	21
ANNEX I: CE CERTIFICATE.....	23

1. USES OF THE INSTRUMENT

Nahita rotary microtome is design with the most advanced technology to slice paraffin included plant or animal tissues in histology routine applications in hospital and research laboratories.


2. DESCRIPTION

Microtome is a very precise instrument used to slice paraffin included tissues. The obtained sections have a micrometer thickness and this fact allows its observation under the microscope.

There are several types of microtomes (hand microtomes, balance, rotary, sliding, criostat, ultramicrotome, etc.) but the most used is the rotary microtome. Some of its advantages are high precision and the possibility of producing very thin serial sections thanks to the desmultiplication process that changes a rotation movement into an advance movement.

The rotary microtome presents the following different parts:

Sample holder: it is the device in which the generally paraffin embedded tissue is placed. Thanks to a pair of screws the paraffin block surface can be oriented in a position parallel to the blade.

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Thank you for choosing this equipment. We sincerely wish that you enjoy your Nahita rotary microtome 591. We highly recommend looking after this equipment according to what is stated in this manual.

Nahita develops its products according to the CE marking regulations as well as emphasizing the ergonomics and security for its user.

The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

The improper use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages, etc. Please read the point of Maintenance, where we expose the security notes.

TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.

Please bear in mind the following:

- ♦ This manual is inseparable from the Nahita rotary microtome 591, so it should be available for all the users of this equipment.
- ♦ You should carefully handle the equipment avoiding sudden movements, knocks, free fall of heavy / sharp objects on it. Avoid spilling liquids inside the equipment.
- ♦ Never dismantle the different pieces of the equipment to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty, as well as injuries on people that handle the microtome.
- ♦ If you have any doubt about setting up, installation or functioning do not hesitate in contacting your wholesaler. You can also tell us any doubts or suggestions you have by contacting Nahita Technical Assistance Department by email to asistencia@auxilab.es or by telephone: +34 807 117 040 (0.30 Euros/min).
- ♦ This equipment is protected under the Warranties and consumer goods regulation (10/2003).
- ♦ Overhaul is not covered by the equipment warranty.
- ♦ Operations made by non-qualified staff will automatically produce a loss of the warranty.
- ♦ Accessories (including their loss) are not covered by the product's warranty. The warranty neither covers piece's deterioration due to the course of time.
- ♦ Please make sure you keep the invoice, either for having the right to claim or asking for warranty coverage. In case you have to send the equipment to Nahita Technical Assistance Department you should enclose the original invoice or a copy as guarantee.
- ♦ Please do not forget filling the warranty certificate and send it before 15 days after the date of purchase.
- ♦ Manufacturer reserves the right to modify or improve the manual or equipment.

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